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Read this manual carefully before use, and retain it for future reference.


Safety Instructions

When using the toothbrush for the first time, it may cause itch or minor gum bleeds. This is normal when the gums are not accustomed to the efficient cleaning of an electric toothbrush. It is recommended to use delicate & gentle mode when your gums are not used to it yet. Normally, gum bleeds will stop within one to two weeks. Stop using the toothbrush and consult a dentist if it persists.

⚠ WARNINGS

- The plug must be removed from the socket-outlet before cleaning or maintaining the appliance.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Do not expose the appliance or battery to excessive temperatures.
- Be aware of the risk of terminals of the battery-operated appliance or battery being short-circuited by metal objects.
- This appliance contains a battery that is non-replaceable.
- Only use and charge the appliance at temperature between 5°C and 35°C.

Product Overview


- | | |
|---------------------------|---------------------|
| 1 Shaft Protective Cover | 2 Shaft |
| 3 Toothbrush Handle | 4 On/Off Button |
| 5 Battery Level Indicator | 6 Mode Indicator |
| 7 Toothbrush Head Cover | • Standard Cleaning |
| 8 Toothbrush Head | • Delicate & Gentle |
| 9 Charging Dock | • Gum Massage |
| 10 User Manual | • Whitening Polish |
| 11 Warranty Notice | |


Note: Illustrations of product, accessories, and user interface in the user manual are for reference purposes only. Actual product and functions may vary due to product enhancements.

How to Use

1. Make sure that the toothbrush head is installed in the same direction as the front of the toothbrush handle, and make sure there is a slight gap between the toothbrush head and the handle (the gap is designed to prevent the toothbrush head from colliding with the handle during high-frequency vibration, so as to generate more power from the motor).

2. Wet the toothbrush head and squeeze a proper amount of toothpaste onto the toothbrush head.

3. Turning on/turning off/changing modes

Turning on: Press the button to turn on the toothbrush.

Changing modes: After the toothbrush is turned on, press within 10 seconds to change modes.


Pause: Press the button for one second to pause the toothbrush. The toothbrush will automatically turn off after being paused for 30 seconds.

Turning off:

- Within 10 seconds after the toothbrush is turned on, press the button to turn it off when the toothbrush is in the whitening polish mode or after changing it to the mode.
- When your toothbrush runs for over 10 seconds, press the button to turn it off under any mode.

4. Healthy brushing

- Position the toothbrush head at an angle of 45° to the gum line and slowly move it along the teeth and gum line.
- During brushing, the toothbrush pauses briefly every 30 seconds to prompt the user to change the brushing position.
- When there is no active pause during use, the toothbrush will automatically turn off after running for 2 minutes.

- Four brushing modes are available:
 - Standard cleaning: Suitable for daily oral hygiene
 - Delicate & gentle: Suitable for new users and people with sensitive teeth
 - Gum massage: Suitable for people with sensitive gums
 - Whitening polish: Suitable for people with stubborn tooth stains

5. Cleaning the toothbrush head

- ➊ After brushing teeth, rinse the toothbrush head with running water.
- ➋ After rinsing the toothbrush head, put on the toothbrush head cover to keep the toothbrush head clean and hygienic.

* Dentists recommend replacing the toothbrush head every three months.


Regular Charging

- Battery level indicator description

Toothbrush Status	Battery Level Indicator Status	Battery Level
Operating	Keeping white	40%–100%
	Keeping orange	20%–40%
	Breathing orange at a slow pace	10%–20%
	Breathing orange at a fast pace (unable to turn on)	<10%
Charging	Keeping white	100%
	Breathing white at a slow pace	40%–100%
	Breathing orange at a slow pace	20%–40%
	Breathing orange at a moderate pace	0–20%

- When the battery level indicator is breathing orange at a fast pace and the toothbrush cannot be turned on, it means the battery is drained. It takes approximately 4 hours (the charging time varies depending on the remaining power) to fully charge the battery. When the battery is drained, the toothbrush can be turned on and complete one brushing session after two minutes of flash charging.

Note: The toothbrush charging dock is designed for this toothbrush and is incompatible with other products.

- Please use a 5 V $\equiv$ 1 A power adapter (sold separately) when charging.

Environmental Protection


When the toothbrush reaches the end of its service life, please remove, dispose of and recycle the built-in battery according to local laws and regulations of the country or region where the product is used. If the battery cannot be removed, dispose of it in strict accordance with the regulations on classified domestic waste disposal in your area.

WARNINGS



- Do not expose the battery to fire to prevent the battery from exploding or releasing toxic and hazardous substances.
- Do not short-circuit the battery to prevent it from combusting and causing a fire.
- Do not disassemble the toothbrush unless you want to dispose of it; otherwise, the toothbrush will be damaged.
- Before disposing of the toothbrush, remove the battery first.
- Turn the toothbrush off before removing the battery.
- The battery should be disposed of safely.
- When the battery is at end of life, the appliance shall be properly disposed of.
- The battery must be removed from the appliance before it is scrapped.
- The appliance must be disconnected from the supply mains when removing the battery.

Removing the rechargeable battery:


1. Insert a flat screwdriver into the gap between the handle's bottom cover and the handle, and pry upward until the bottom cover opens.

2. Use a flat screwdriver to press the buckles on both sides of the inner compartment, press the front metal shaft against a hard object, and push the inner compartment out of the toothbrush.

3. Use a diagonal cutter to cut off the positive and negative battery connection wires. Remove the battery and deposit it at an official recycling point.

Specifications

Name: Xiaomi Electric Toothbrush T302	
Model: MES608	
Rated Voltage: 3.6 V $\overline{\sim}$	Rated Power: 2 W
IP Rating: IPX8	Battery: 3.6 V $\overline{\sim}$, 1050 mAh, 3.78 Wh
Charging Time: Approx. 4 hours	

Name: Charging Dock	
Model: MES608	Rated Input: 5 V $\overline{\sim}$ 1 A
Charging Method: Wireless charging	
Maximum Output Power: -13 dBm	
Operating Frequency: 110–205 kHz	
This device has been tested and meets applicable limits for radio frequency (RF) exposure.	

Troubleshooting

Issues	Possible causes	Solutions
The toothbrush does not turn on.	The toothbrush cannot be turned on while charging.	Unplug the charging cable from the toothbrush.
	The battery level is too low.	Charge the battery in a timely manner.
The toothbrush head cannot be installed.	There are foreign objects inside the toothbrush head.	Use a tool to remove the foreign objects from the toothbrush head.
	There are foreign objects on the shaft.	Clean the shaft.
	The toothbrush head is not a genuine product.	Replace it with a genuine original toothbrush head.
The toothbrush head vibration is weak.	The toothbrush head is not installed in place.	Reinstall the toothbrush head.
	There are foreign objects between the toothbrush head and the handle.	Remove the toothbrush head and remove the foreign objects.

WEEE Disposal and Recycling Information

All products bearing this symbol are waste electrical and electronic equipment (WEEE as in directive 2012/19/EU) which should not be mixed with unsorted household waste. Instead, you should protect human health and the environment by handing over your waste equipment to a designated collection point for the recycling of waste electrical and electronic equipment, appointed by the government or local authorities. Correct disposal and recycling will help prevent potential negative consequences to the environment and human health. Please contact the installer or local authorities for more information about the location as well as terms and conditions of such collection points.

Manufactured for: Xiaomi Communications Co., Ltd.
Manufactured by: Shenzhen Soocas Technology Co., Ltd.
(a Mi Ecosystem company)

Address: Room 401, 402, 1101, 1102, 1103, 1104, Building 2 Chongwen Area, Nanshan iPark, No. 3370 Liuxian Avenue, Fuguang Community, Taoyuan Subdistrict, Nanshan District, Shenzhen, China
For further information, please go to www.mi.com
Made in China

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Fabricado para: Xiaomi Communications Co., Ltd.
Fabricado por: Shenzhen Soocas Technology Co., Ltd.